

TED HOSE SIZING GUIDE


- 1. Choose between Knee High or Thigh High Ted Hose
- 2. Measurements needed for proper sizing:


Calf Circumference

Measure at widest portion of calf

Length

Knee High measure heel to fold of knee
*Thigh High measure heel to gluteal fold





Sizing Chart

MSC			Calf	
Order #	<u>Style</u>	Size	Circumference	Length
2501	Knee	Small Regular	up to 12"	up to 16"
2502	Knee	Small Long	up to 12"	greater than 16"
2503	Knee	Medium Regular	between 12-15"	up to 17"
2504	Knee	Medium Long	between 12-15"	greater than 17"
1748	Knee	Large Regular	between 15-17 1/2"	up to 18"
1749	Knee	- Large Long	between 15-17 1/2"	greater than 18"
20091	Knee	X-Large Regular	greater than 17 1/2"	up to 18"
20090	Knee	X-Large Long	greater than 17 1/2"	greater than 18"
MSC			Calf	
Order #	<u>Style</u>	Size	Circumference	Length
20097	Thigh	Small Short	up to 12"	up to 29"
20094	Thigh	Small Regular	up to 12"	between 29-33"
20093	Thigh	Small Long	up to 12"	greater than 33"
20096	Thigh	Medium Short	between 12-15"	up to 29"
20100	Thigh	Medium Regular	between 12-15"	between 29-33"
20092	Thigh	Medium Long	between 12-15"	greater than 33"
20095	Thigh	Large Short	greater than 15"	up to 29"
20099	Thigh	Large Regular	greater than 15"	between 29-33"
20098	Thigh	Large Long	greater than 15"	greater than 33"
LIMOPC 2850 South Industrial Highway - Suite 400				

UMOPC 2850 South Industrial Highway • Suite 400 • Ann Arbor, Michigan 48104 • (734) 973-2400 University Hospital Satellite • 1500 E. Medical Center Drive • Ann Arbor, Michigan 48109 • (734) 936-7043